

PTA Balloon Catheter


0.035"

The PTA Balloon Catheter is a **high performance** balloon catheter for peripheral indication. The device features an **ultra-low profile, semi-compliant** balloon combined with a low profile tip.

Our focus is to enhance technology by developing devices with Low Entry Tip Profiles, guided by Ergonomic Pushability for your ease of use, especially in complex vessels, Tactile Trackability for enhanced Manoeuvring.

As an essential feature the Dual Lumen Shaft Technology provides a rapid inflation/deflation time.

	Competitor A	Competitor B	Competitor C	PTA 0.035"
TIP ENTRANCE DIAMETER compatibility comparison	0.048"	0.040"	0.048"	0.037"




PTA Balloon Catheter

0.035"

TECHNICAL SPECIFICATIONS

Description	PTA Balloon Catheter Semi-Compliant OTW
Recommended Guidewire	0.035" (0.89 mm)
Recommended Introducer Sheath	5 F - 9 F
Nominal Pressure	3 mm to 12 mm - 7 bar
Rated Burst Pressure (RBP)	<ul style="list-style-type: none"> • Ø 3 mm - 13 bar • Ø 3.5 mm to 12 mm - 10 bar
Working Catheter Length	80 cm 120 cm 150 cm
Platinum X-ray Marker	One marker at proximal and distal end of balloon
Low Entry Profile	0.95 mm (0.037")
Balloon Material	Pebax
Hydrophilic Coating	400 mm distal

ORDER INFORMATION

Product Code 80 cm length	Product Code 120 cm length	Product Code 150 cm length	Balloon Diameter (mm)	Balloon Length (mm)	Minimum Balloon Profile (mm)	Recommended Introducer Sheath (F)
08PVQ100010APO	12PVQ100010APO	15PVQ100010APO	10	10	1.79	6
08PVQ100015APO	12PVQ100015APO	15PVQ100015APO	10	15	1.79	6
08PVQ100020APO	12PVQ100020APO	15PVQ100020APO	10	20	1.79	6
08PVQ100030APO	12PVQ100030APO	15PVQ100030APO	10	30	1.79	6
08PVQ100040APO	12PVQ100040APO	15PVQ100040APO	10	40	1.79	6
08PVQ100060APO	12PVQ100060APO	15PVQ100060APO	10	60	1.79	6
08PVQ100080APO	12PVQ100080APO	15PVQ100080APO	10	80	1.79	6
08PVQ100100APO	12PVQ100100APO	15PVQ100100APO	10	100	1.79	6
08PVQ100120APO	12PVQ100120APO	15PVQ100120APO	10	120	1.79	6
08PVQ100150APO	12PVQ100150APO	15PVQ100150APO	10	150	1.79	6
08PVQ100200APO	12PVQ100200APO	15PVQ100200APO	10	200	1.79	6
08PVQ100220APO	12PVQ100220APO	15PVQ100220APO	10	220	1.79	6
08PVQ100280APO	12PVQ100280APO	15PVQ100280APO	10	280	1.79	6
08PVQ120010APO	12PVQ120010APO	15PVQ120010APO	12	10	2.16	9
08PVQ120015APO	12PVQ120015APO	15PVQ120015APO	12	15	2.16	9
08PVQ120020APO	12PVQ120020APO	15PVQ120020APO	12	20	2.16	9
08PVQ120030APO	12PVQ120030APO	15PVQ120030APO	12	30	2.16	9
08PVQ120040APO	12PVQ120040APO	15PVQ120040APO	12	40	2.16	9
08PVQ120060APO	12PVQ120060APO	15PVQ120060APO	12	60	2.16	9
08PVQ120080APO	12PVQ120080APO	15PVQ120080APO	12	80	2.16	9
08PVQ120100APO	12PVQ120100APO	15PVQ120100APO	12	100	2.16	9
08PVQ120120APO	12PVQ120120APO	15PVQ120120APO	12	120	2.16	9
08PVQ120150APO	12PVQ120150APO	15PVQ120150APO	12	150	2.16	9
08PVQ120200APO	12PVQ120200APO	15PVQ120200APO	12	200	2.16	9
08PVQ120220APO	12PVQ120220APO	15PVQ120220APO	12	220	2.16	9
08PVQ120280APO	12PVQ120280APO	15PVQ120280APO	12	280	2.16	9

CE 0124



QualiMed Innovative Medizinprodukte GmbH
 Boschstraße 16, 21423 Winsen, Germany
 Phone +49 4171 6578 0
 Fax +49 4171 6578 11
 Web www.qualimed.de

Q Medical GROUP
DEVICES LIMITED


BRO_PTA035 Rev A

Peripheral Vascular



PTA Balloon Catheter

0.035"



Low-Profile

Regulatory Status
 CE Mark Approved



